

## Solvent Data Chart

# NMR

Acetic Acid-d<sub>4</sub>  
Acetone-d<sub>6</sub>  
Acetonitrile-d<sub>3</sub>  
Benzene-d<sub>6</sub>  
Chloroform-d  
Cyclohexane-d<sub>12</sub>  
Deuterium Oxide  
N, N-Dimethylformamide-d<sub>7</sub>  
Dimethyl Sulfoxide-d<sub>6</sub>  
1,4-Dioxane-d<sub>8</sub>  
Ethanol-d<sub>6</sub>  
Methanol-d<sub>4</sub>  
Methylene Chloride-d<sub>2</sub>  
Pyridine-d<sub>5</sub>  
1,1,2,2-Tetrachloroethane-d<sub>2</sub>  
Tetrahydrofuran-d<sub>8</sub>  
Toluene-d<sub>8</sub>  
Trifluoroacetic Acid-d  
Trifluoroethanol-d<sub>3</sub>

<sup>1</sup> H Chemical Shift (ppm from TMS) (multiplicity)	JHD (Hz)	<sup>13</sup> C Chemical Shift (ppm from TMS) (multiplicity)	JCD (Hz)	<sup>1</sup> H Chemical Shift of HOD (ppm from TMS)	Density at 20°C	Melting point (°C)	Boiling point (°C)	Dielectric Constant	Molecular Weight
11.65 (1) 2.04 (5)	2.2	178.99 (1) 20.0 (7)	20	11.5	1.12	16.7	118	6.1	64.08
2.05 (5)	2.2	206.68(1) 29.92 (7)	0.9 19.4	2.8*	0.87	-94	56.5	20.7	64.12
1.94 (5)	2.5	118.69 (1) 1.39(7)	21	2.1*	0.84	-45	81.6	37.5	44.07
7.16 (1)		128.39(3)	24.3	0.4	0.95	5.5	80.1	2.3	84.15
7.24 (1)		77.23	32.0	1.5*	1.50	-63.5	61-62	4.8	120.38
1.38 (1)		26.43 (5)	19	0.8	0.89	6.47	80.7	2.0	96.24
4.80 (DSS) 4.81 (TSP)		NA	NA	4.8	1.11	3.81	101.42	78.5	20.03
8.03 (1) 2.92 (5) 2.75 (5)	1.9 1.9	163.15 (3) 34.89 (7) 29.76 (7)	29.4 21.0 21.1	3.5	1.03	-61	153	36.7	80.14
2.50 (5)	1.9	39.51 (7)	21.0	3.3*	1.19	18.55	189	46.7	84.17
3.53 (m)		66.66 (5)	21.9	2.4	1.13	11.8	101.1	2.2	96.16
5.19 (1) 3.56 (1) 1.11(m)		56.96 (5) 17.31 (7)	22 19	5.3	0.89	-114.1	78.5	24.5	52.11
4.78 (1) 3.31 (5)	1.7	49.15 (7)	21.4	4.9	0.89	-97.8	64.7	32.7	36.07
5.32 (3)	1.1	54.00 (5)	27.2	1.5	1.35	-95	39.75	8.9	86.95
8.74 (1) 7.58 (1) 7.22 (1)		150.35 (3) 135.91 (3) 123.87 (3)	27.5 24.5 25	5	1.05	-41.6	115.2-115.3	12.4	84.13
6.0		73.78 (3)			1.62	-44	146.5	8.20	169.86
3.58 (1) 1.73 (1)		67.57 (5) 25.37 (5)	22.2 20.2	2.4-2.5	0.99	-108.5	66	7.6	80.16
7.09 (m) 7.00 (1) 6.98 (5) 2.09 (5)	2.3	137.86 (1) 129.24 (3) 128.33 (3) 125.49 (3) 20.4 (7)	23 24 24 19	0.4	0.94	-95	110.6	2.4	100.19
11.50 (1)		164.2 (4) 116.6 (4)		11.5	1.49	-15.4	72.4		115.03
5.02 (1) 3.88 (4x3)	2(9)	126.3 (4) 61.5 (4x5)	22	5	1.41	-43.5	74.05		103.06

M.J. O'Neil, P.E. Heckelman, C.B.Koch, K.J. Roman, *The Merck Index*, an Encyclopedia of Chemicals, Drugs, and Biologicals - Fourteenth Edition, Merck Co., Inc. Whitehouse Station, NJ 2006.

- The <sup>1</sup>H spectra of the residual protons and <sup>13</sup>C spectra were obtained on a Varian Gemini 200 spectrometer at 29.5°K. The NMR solvents used to acquire these spectra contain a maximum of 0.05% and 1.0% TMS (v/v) respectively. Since deuterium has a spin of 1, triplets arising from coupling to deuterium have the intensity ratio of 1:1:1. 'm' denotes a broad peak with some fine structures. It should be noted that chemical shifts can be dependent on solvent, concentration and temperature.
- Approximate values only, may vary with pH, concentration and temperature.
- Melting and boiling points are those of the corresponding unlabeled compound (except for D<sub>2</sub>O). These temperature limits can be used as a guide to determine the useful liquid range of the solvents. Information gathered from the Merck Index -14th Edition.
- ★ HOD Peaks - NMR spectra of "neat" deuterated solvent always exhibit a peak due to H<sub>2</sub>O in addition to the residual solvent peak. When the exchange rate between H<sub>2</sub>O and HDO is slow on the NMR timescale the water peak appears as two peaks, a singlet corresponding to H<sub>2</sub>O and a 1:1:1 triplet corresponding to HDO.