Pesticide Residue on Organic Fruit Vs. Non-Organic Fruit Alexandra Heuston, B.S. Candidate in Chemistry California State University Stanislaus

Introduction

- Pesticides are necessary in agriculture to protect crops from Insect infestation, diseases and weed takeover.
- The United States has maximum residue levels (MRLs) set by the EPA which are enforced and monitored by the USDA.
- Americans trust the USDA to ensure the food we consume daily is within a healthy standard.
- Fruits and vegetables are categorized as organic and conventional/ non-organic.
- It is a common myth that organic foods do not contain pesticides, a myth most consumers believe due to miseducation

Background & Literature Review

- Montiel et. al analyzed organic and conventional produce using QuEChERS with LC-MS/MS in Quebec, Canada and found higher level of detected pesticides for conventional than organics produce (Montiel, 2019).
- Tobin *et. al* similarly found in Ireland one or more pesticide residues in 15 out of 27 organic samples. In conventional samples, 17 had pesticide residues present out of 27 (Tobin, 2014).
- Bursi'c *et al.* in Siberia also tested an identical question sampling different produce from around Europe and found high pesticide residue for both and found 9 MRL violations on both organic and conventional apples (Bursi'c, 2021).
- An independent study has not been done in the United States outside of the USDA who
 is not region specific with their sampling.
- The USDA tested organic fruits and vegetables in 2010 and found 21 samples out of 571 in violation of USDA organic regulations and none in violation of EPA tolerances (USDA, 2012).
- The USDA uses LC-MS/MS for detection, identification and quantification. DART, a new method developed in 2005, proves by Farré *et al.* to be a novel method for quick food quality detection (Farré, 2013).

Research Questions

What identity of pesticides can be found on organic fruit in comparison to conventional fruit within the Central Valley of California? What is the concentration of each pesticide found on various fruit samples in comparison to EPA tolerances?

Method

A quantitative research approach will be implemented to detect and identify pesticides remaining on fruit and to quantify the concentrations of the pesticides found to compare against EPA MRLs.

Samples

- Strawberries and apples are the commodities of choice. Strawberries are the #1 fruit with pesticide violations and apples are #5 (USDA, as cited by EWG 2021). Strawberries and apples are the 2nd and 3rd most bought fruits respectively after bananas (The Packer Fresh Trends, 2019).
- A standard container of USDA organic and conventional strawberries as well as a presorted bundle of USDA organic apples and conventional apples purchased at a local grocer in Stanislaus County will be compared.

Instrumentation

- DART-MS/MS coupled to a linear quadrupole ion trap will be used for direct, rapid detection of pesticides due to little same preparation required, to screen for pesticides presence
- The USDA method of detection and quantification will then be used after DART. The QuEChERS method will be used to extract and isolate the pesticides residing on the samples, adapted from (Anastassiades, 2003). HPLC (High-Performance Liquid Chromatography) will then be used to perform analysis.

Expected Results

- The number of pesticides detected residually on organic fruit will be less than the number found for conventional fruits.
- The concentration of each pesticide found for organic fruits will be less than the concentrations found on conventional fruit.
- Strawberries will have more pesticides detected than apples

Significance

- The research aims to inform this region, the Central Valley of California, on the certain pesticides that reside on the fruit they consume.
- The public will be informed on how the levels of detected pesticides lie in terms of concentration compared to EPA MRLs.

References

- Anastassiades M, Lehotay SJ, Stajnbaher D, Schenck FJ. Fast and easy multiresidue method employing acetonitrile extraction/partitioning and "dispersive solidphase extraction" for the determination of pesticide residues in produce. J AOAC Int. 2003 Mar-Apr;86(2):412-31. PMID:12723926
- 2. Bursi'c, V.; Vukovi'c, G.; Cara, M.; Kosti'c, M.; Stojanovi'c, T.; Petrovi'c, A.; Puva'ca, N.; Marinkovi'c, D.; Konstantinovi'c, B. Plant Protection Products Residues
 Assessment in the Organic and Conventional Agricultural Production.
 Sustainability 2021, 13, 1075. https://doi.org/10.3390/su13031075
- Farré, M., Picó, Y., & Barceló, D. (2013). Direct Peel Monitoring of Xenobiotics in Fruit by Direct Analysis in Real Time Coupled to a Linear Quadrupole Ion Trap—Orbitrap Mass Spectrometer. *Analytical Chemistry*, 85(5), 2638–2644. https://doi.org/10.1021/ac3026702
- 4. Farm Journal Media. (2019, March 28). The Packer Fresh Trends 2019. http://digitaledition.qwinc.com/publication/?m=40749&i=57744&p=1.

- Group, E. W. (n.d.). EWG's 2021 Shopper's Guide to Pesticides in Produce [™]. EWG's 2021 Shopper's Guide to Pesticides in Produce | Summary. https://www.ewg.org/foodnews/summary.php.
- Montiel-León, J. M., Duy, S. V., Munoz, G., Verner, M.-A., Hendawi, M. Y., Moya, H., ...
 Sauvé, S. (2019, April 26). Occurrence of pesticides in fruits and vegetables from
 organic and conventional agriculture by QuEChERS extraction liquid
 chromatography tandem mass spectrometry. Food Control.
 https://www.sciencedirect.com/science/article/abs/pii/S0956713519301884
- 7. Tobin R, Walsh T, Garvey J, Larkin T (2014) Detection of Pesticide Residues in Organic and Conventional Fruits and Vegetables Available in Ireland Using Gas Chromatography/Tandem Mass Spectrometry (GC-MS/MS) and Liquid Chromatography/Tandem Mass Spectrometry (LC-MS/MS) Detection. J Nutrition Health Food Sci 2(1): http://dx.doi.org/10.15226/jnhfs.2014.00
- 8. USDA National Organic Program. (2012, November 1). 2010-2011 Pilot Study: Pesticide Residue Testing of Organic Produce. USDA Agricultural Marketing Services. https://www.ams.usda.gov/reports/2010-2011-pilot-study-pesticide-residue-testing-organic-produce